pós-graduação em ciências da saúde do iamspe


Genômica Computacional

Analise de variantes

Professor: Ricardo A. Vialle

CS31 - Genômica Computacional


Read mapping



Read mapping

```
221
         231
                   241
                            251
                                     261
                                             271
                                                       281
                                                                291
                                                                          301
                                                                                   311
                                                                                            321
TCGAGC*GCCATCTGCGGGCCC
                                 tcagccgccgactcggacatcagatccacctcggggtccgc*atcgctccgccggcgtccgacgaagccgacgac
TCGAGC*GCCATCTG
                                     ccgccgactcggacatcagatccaactcggggtccgc*atcgctccgccggcgaccgacgaagccgacga
TCGAGC*GCCAT
                                             cggacatcagatccaactcggggtccgc*atcgctccgtcggcgaccgacgaagccgacga
TCGAGC*GCCATCTGCGGGC
                                             cggacatcagatccaactcggggtccgc*atcgctccgccggcgaccgacgaagccgacga
tca
           CTGCGGGCCCCGATTACCGACGCTCTCAGCCGCCGACT
                                             CGGACATCAGATCCAACTCGGGGTCCGC*ATCGCTCCGCCGGCGACCGACGAAGCCGACGAI
tcgagc*gccatctgcgggccc
                                               GACATCAGATCCAACTCGGGGTCCGC*ATCGCTCCGCCGGCGACCGACGAAGCCGACGA
tcgagc*gccatctgcgggccccgattaccgacgctctcag
                                                CATCAGATCCAACTCGGGGTCCGC*ATCGCTCCGCCGGCGACCGACGAAGCCGACGA
TCGAGC*GCCATCTGCGGGCCCCGATTACCGACGCTCTCAGCCGCCGAC
                                             CGG
                                                  tcagatccaactcggggtccgc*atcgctccgccggcgaccgacgaagccgacga
                                             CGG
                                                         AACTCGGGGTCTGC*ATCGCTCTGCCGGCGACCGACGAGCCGACGA
TCGAGC*GCCATCTGCGGGCCCCGATTACCGACGCTCTCAGCCGCCGAC
TCGAGC*GCCATCTGCGGGCCCCGATTACCGAAGCTCTCAGCCGCCGAC
                                            CGGA
                                                             cgggggccgc*atcgcgccgtcggcgaccgacgaagccgacac
                                                                            CCGCCGGCGACGACGAAGCCGACGAI
TCGAGC*GCCATCTGCGGGCCCCGATTACCGACGCTCTCAGCCGCCGAC
                                             CGGACAT
TCGAGC*GCCATCTGCGGGCCCCGATTACCGACGCTCTCAGCCGCCGACTCGGACATCAG
                                                                              cccggcgaccgacaaagccgacga
TCGAGC*GCCATCTGCGGGCCCCGATTACCGACGCTCTCAGCCGCCGACTCGGACATCAG
                                                                                             GCCGACGA
tcgagc*gccatctgcgggccccgattaccgacgctctcagccgccgactcggacatcaga
                                                                                             gccgacga
tcgagc*gccatctgcgggccccgattaccgacgctctcagccgccgactcggacatcaga
                                                                                             ccgacga
TCGAGC*GCCATCTGCGGGCCCCGATTACCGACGCTCTCAGCCGCCGACTCGGACATCAGATCCA
                                                                                                acga
TCGAGC*GCCATCTGCGGGCC
                                                                                                 CGA
tcgagc*gccatctgcgggcc
tcgagc*gccatctgcgggccccgattaccgacgctctcagccgccgactcggacatcagatc
TCGAGC*GCCATCTGCGCGCCCCGATTACCGACCC
TCGAGC*GCCATCTGCGGGCCCCGATTACCGACGCTCTCAGCCGCCGACTCGGACATCAGATCCAACTCGGGGTCCGC*AT
TCGAGC*GCCATCTGCGGGCCCCGATTACCGACGCTCTCAGCCGCCGACTCGGACATCAGATCCAACTCGGGGTCCGC*ATCGCTCCGCCGGCGACCGA
TCGAGC*GCCATCTGCGGGCCCCGATTACCGACGCTCTCAGCCGCCGACTCGGACATCAGATCCAACTCGGGGTCCGC*A
```

IGV browser



Types of genetic variation


- 99% of DNA is shared between two individuals
- Variation in the remainder explains all our predisposition differences
- Remaining phenotypic variation: environmental/stochastic differences

Name	Example	Frequency in one genome
Single nucleotide polymorphisms (SNPs)	GAGGAGAACG[C/G]AACTCCGCCG	1 per 1,000 bp
Insertions/deletions (indels)	CACTATTC[C/CTATGG]TGTCTAA	1 per 10,000 bp
Short tandem repeats (STRs)	ACGGCAGTCGTCGTCACCGTAT	1 per 10,000 bp
Structural variants (SVs) / Copy Number Variants (CNVs)	Large (median 5,000 bp) deletions, duplications, inversions	1 per 1,000,000 bp

VCF and BCF file format

VCF

```
##fileformat=VCFv4.2
 ##contig=<ID=2,length=51304566>
 ##INFO=<ID=AC, Number=A, Type=Integer, Description="Allele count in genotypes">
 ##INFO=<ID=AN, Number=1, Type=Integer, Description="Total number of alleles in called genotypes">
 ##FORMAT=<ID=GT, Number=1, Type=String, Description="Genotype">
 ##FORMAT=<ID=DP, Number=1, Type=Integer, Description="Read Depth">
 ##FORMAT=<ID=GQ, Number=1, Type=Integer, Description="Genotype Quality">
 #CHROM POS ID REF ALT QUAL FILTER INFO FORMAT
                                                      SAMPLE1
                                                                  SAMPLE2
                                                                              SAMPLE3
                                                                                          SAMPLE4
                                                                                                      SAMPLE5
                                                                                                                  SAMPLE6
                                                                                                                               SAMPLE7
 2 81170
          . C T
                            AC=9; AN=7424
                                           GT:DP:GQ
                                                      0/0:4:12
                                                                  0/0:3:9
                                                                              0/1:1:3
                                                                                          0/1:9:24
                                                                                                      1/0:4:12
                                                                                                                  0/0:5:15
                                                                                                                               0/0:4:12
          . G A
                                                                              0/0:1:3
                                                                                          0/0:9:24
                                                                                                      0/1:4:12
                                                                                                                               0/0:4:12
   81171
                            AC=6; AN=7446
                                           GT:DP:GQ
                                                      0/1:4:12
                                                                  0/0:3:9
                                                                                                                  0/1:5:15
   81182
             A G
                            AC=5; AN=7506
                                           GT:DP:GQ
                                                      0/0:5:15
                                                                  0/0:4:12
                                                                              0/0:5:15
                                                                                          0/0:9:24
                                                                                                      0/0:4:12
                                                                                                                  0/0:4:12
                                                                                                                               0/0:4:12
 2 81204
           . T G
                                                                              0/0:10:30
                                                                                          0/0:15:39
                                                                                                      0/0:9:27
                                                                                                                  1/0:13:39
                                                                                                                               0/1:14:42
                            AC=2; AN=7542
                                           GT:DP:GQ
                                                      1/0:5:15
                                                                  0/0:9:27
BCF
                                                                                                       GQ:12: 9: 3:24:12:15:12
           . C T
                                                                             DP:4:3:1:9:4:5:4
   81170
                            AC=9; AN=7424
                                           GT:0/0:0/0:0/1:0/1:1/0:0/0:0/0
          . G A
    81171
                            AC=6; AN=7446
                                           GT:0/1:0/0:0/0:0/0:0/1:0/1:0/0
                                                                             DP:4:3:1:9:4:5:4
                                                                                                       GQ:12: 9: 3:24:12:15:12
    81182
          . A G
                            AC=5; AN=7506
                                           GT:0/0:0/0:0/0:0/0:0/0:0/0:0/0
                                                                             DP:5:4:5:9:4:4:4
                                                                                                       GQ:15:12:15:24:12:12:12
 2 81204
          . T G
                            AC=2; AN=7542
                                           GT:1/0:0/0:0/0:0/0:0/0:1/0:0/1
                                                                             DP:5:9:10:15:9:13:14
                                                                                                       GQ:15:27:30:39:27:39:42
```




QNAME: query template name, aka. read ID

SAM and BAM file format (FLAG)

Bitwise Flags

Integer	Binary	Description (Paired Read Interpretation)
1	000000000001	template having multiple templates in sequencing (read is paired)
2	000000000010	each segment properly aligned according to the aligner (read mapped in proper pair)
4	000000000100	segment unmapped (read1 unmapped)
8	000000001000	next segment in the template unmapped (read2 unmapped)
16	000000010000	SEQ being reverse complemented (read1 reverse complemented)
32	000000100000	SEQ of the next segment in the template being reverse complemented (read2 reverse complemented)
64	000001000000	the first segment in the template (is read1)
128	000010000000	the last segment in the template (is read2)
256	000100000000	not primary alignment
512	001000000000	alignment fails quality checks
1024	010000000000	PCR or optical duplicate
2048	100000000000	supplementary alignment (e.g. aligner specific, could be a portion of a split read or a tied region)

Use *samtools flagstat* to summarize these metrics



QNAME: query template name, aka. read ID

CIGAR

For example:

```
RefPos: 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 Reference: C C A T A C T G A A C T G A C T A A C
```

Read: ACTAGAATGGCT


Aligning these two:

```
RefPos: 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 Reference: C C A T A C T G A C T A A C Read: A C T A G A A T G G C T
```


With the alignment above, you get:

POS: 5


CIGAR: 3M1I3M1D5M



QNAME: query template name, aka. read ID



The smallest template-independent mapping quality in the template Alignment score generated by aligner Barcode sequence identifying the sample Offset to base alignment quality (BAQ) Phred quality of the unique molecular barcode bases in the OX tag	SA:Z:ref,29,-,6
Barcode sequence identifying the sample Offset to base alignment quality (BAQ)	
Offset to base alignment quality (BAQ)	
Phred quality of the unique molecular barcode bases in the OX tag	= SA · Z · ref 9 + 55
Cell identifier	- NM:i:1
Reference name of the next hit	
BAM only: CIGAR in BAM's binary encoding if (and only if) it consists of >65535 operators	Optional fields in
Edit distance between the color sequence and the color reference (see also NM)	EUAL: read quality; * me
Free-text comments	CHI
Leftmost coordinate of the next hit	leftmost/rightmost read.
Color read base qualities	griment of the NEXT rea
	BAM only: CIGAR in BAM's binary encoding if (and only if) it consists of >65535 operators Edit distance between the color sequence and the color reference (see also NM) Free-text comments Leftmost coordinate of the next hit



RNEXT: reference sequence name of the primary alignment of the NEXT read. For paired-end sequencing, NEXT read is the paired read, corresponding to the RNAME column.

CIGAR: summary of alignment, e.g. insertion, deletion

MAPQ: mapping quality

POS: 1-based position


RNAME: reference sequence name, e.g. chromosome/transcript id

FLAG: indicates alignment information about the read, e.g. paired, aligned, etc.


QNAME: query template name, ake, read ID



PRE-PROCESSING



VARIANT DISCOVERY



CALLSET REFINEMENT

